

ELECTRONIC DENSITY (SPECIFIC GRAVITY) BALANCE



DETAILS:

This is a special balance to measure solid (liquid) density (specific gravity), which is widely applied to measure resin, metal, glass and chemical liquid in electronics, rubber, food, paper and machinery processing industries.

Specific gravity is measured in immersion method under Archimedes' principle. Operate as per clues displayed on the balance, then specific gravity value can be directly, rapidly and conveniently read out.

MODEL:

PLBA 1097

PLBA 1098

	PLBA 1097	PLBA 1098
Max Range	100G	200G
Min Range	10Mg	
Fixed Support	H=203MM	
Measuring Platform	Φ110mm	
Output	RS-232C	
Net Weight	1Kg	
Product Size	250*180*246mm	
Package Size	335*255*252mm	
Accessories	Measuring platform, measuring cup, fixed support, thermometer, doublehook weights, platinum wire, plumb bob, wire basket, balance body, adaptor and standard weights.	