

Medical Refrigerator

Medical Refrigerator is a specialized refrigerating equipment for cold storage of pharmaceuticals as well as biological products and so forth. It is suitable for hospitals, drugstore, pharmaceutical factories, sanitation and antiepidemic stations and clinics.



50L



120L&250L&300L



450L



650L



1000L



1400L

Technical Parameter:

Model		PLC050	PLC051	PLC052		PLC053	PLC054	PLC055	PLC056	PLC057
External Size(L*W*H)		430*450*575mm	450*430*1460mm	580*	560*1785mm	620*560*1970mm	660*700*1980mm	1100*650*2050mm	1260*700*1980mm	1860*760*1980mm
Capacity		50L	120L	250L		300L	450L	650L	1000L	1400L
Temperature Range		2℃ ~ -8℃								
Temperature Accuracy		0.1℃						±2℃	0.1℃	
Control System		Microprocessor control								
Display		Large screen LED display								
Alarm		Audio and visual alarm for: High and low temperature; System & Sensor failure; Door ajar(except PLC050) ; Power failure alarm								
Refrigeration Type		Forced air refrigeration system								
Refrigerant		R134a, CFC-Free								
Compressor		International famous compressor								
Condenser & Evaporator		High efficiency air-cooled condenser, made of cupreous.								
Defrost		Auto defrost								
Structure		Unibody Design & Mono-assembly foaming (Rigid polyurethane insulation material)								
Material	Internal	ABS plate	Cold-rolled steel with anti-bacteria powder coating							
	External	Cold-rolled steel with anti-bacteria powder coating								
Door		Foam	Three-ply transparent	glass door with lock; Heating film preventing frosting; Sealing joint strip with magnetic force						
Shelves		Adjustable shelves made of quality steel wire								
		1	3	4	4	5	10	16	24	
Casters		No	4 universal casters with brake							
Consumption		50W	93W	176W		193W	450W	550W	550W	650W
Power Supply		110/220V±10%, 50/60Hz						175/215		
Standard Accessory		Shelves	Fluorescent lamp, Shelves					shipping dimation: 1150*720*2150		