

## Horizontal Laminar Flow Cabinet



<b>Model</b>	PLS359	PLS360
<b>External Size(W*D*H)</b>	1060*800*1690mm	1460*800*1690mm
<b>Internal Size(W*D*H)</b>	980*500*600mm	1380*500*600mm
<b>Work Surface Height</b>	700mm	
<b>Display</b>	LED Display	
<b>Airflow Velocity</b>	Average of 0.3~0.5m/s	
<b>Material</b>	Main Body:Cold-rolled steel with anti-bacteria powder coating	
	Work Table:304 stainless steel	
	Side Window:5mm toughened glass,anti-ultraviolet radiation	
<b>Pre-filter</b>	Polyester fiber, washable	
<b>HEPA Filter</b>	99.999% efficiency at 0.3um	
<b>Noise</b>	<60 dB	
<b>Fluorescent Lamp</b>	21W*1	28W*1
<b>UV Lamp</b>	20W*1	30W*1
<b>Consumption</b>	320W	620W
<b>Caster</b>	Universal wheel with leveling feet	
<b>Power Supply</b>	110/220V±10%;50/60Hz	
<b>Standard Accessory</b>	Fluorescent lamp,UV lamp*2,Base stand,Waterproof socket*2	
<b>Gross Weight</b>	316 lbs	404 lbs
<b>Package Size(W*D*H)</b>	1230*940*1900mm	1630*940*1900mm