

## Biochemistry Incubator

### Features:

1. Stainless steel inner chamber, round angle structure, adjustable shelves.
2. Microprocessor control system, LCD display.
3. CFC-free refrigeration and gentle breeze circulation.
4. Double-layer glass observation window.
5. Function of self-diagnosis, sensor fault alarm and over-temperature protection.

### Optional accessories

- a. Desk printer
- b. RS485 interface for software connection
- c. BOD sockets
- d.  $\phi$ 50mm test port.



### Technical Parameters:

Model	PLT100	PLT101	PLT102	PLT103	PLT104
Capacity	70L	100L	150L	250L	150L
Temperature Range	0~60°C				-8~60°C
Temperature Precision	0.1°C				
Temperature Fluctuation	±0.5°C				
Ambient Temperature	+4~35°C				
Consumption	350W	400W	500W	680W	650W
Power Supply	110/220V±10%,50/60Hz				
Timing Range	1~9999mins				
Shelves No.	2pcs	2pcs	3pcs	3pcs	3pcs
Internal Size(W*D*H)(mm)	450*330*450	510*383*505	500*370*800	550*445*1100	500*370*800
External Size(W*D*H)(mm)	560*535*798	610*585*930	660*620*1465	650*695*1752	600*620*1465
Package Size(W*D*H)(mm)	640*620*950	690*670*1080	740*700*1710	760*850*1920	740*700*1710
Gross Weight	187.5 lbs	231.5 lbs	291 lbs	326.5 lbs	291 lbs