

LIQUID DENSITY (SPECIFIC GRAVITY) BALANCE



DETAILS:

This balance is widely applied to measure resin, metal, glass and chemical liquid in electronics, rubber, food, paper, and machinery processing industries.

The adoption of agate cushion guarantee good flexibility and accuracy of the balance.

Substantially package is used for good portability and easy operation.

Key parts are made of stainless steel with high intensity and good stability.

MODEL: **PLBA 1096**

Specific gravity range	0~2.000
Standard	1.00 for the pure water in temperature of 20°C
Plumb bob discharge capacity	5cm ³
Accuracy	0.001
External size	215*90*262 mm
Package size	255*186*133 mm
Net weight	3.3 Lbs