

Technical Parameters:

Model	PLS200	PLS201	PLS202	PLS203	PLS204
Internal Size (W*D*H)(mm)	940*600*660	1100*600*660	1300*700*630	1710*600*660	1800*600*660
External Size (W*D*H)(mm)	1100*750*2250	1300*850*2200	1500*840*2280	1870*750*2290	1950*800*2305
Tested Opening	Safety Height \leq 200mm				
Max Opening	400mm				
Inflow Velocity	0.53 \pm 0.025m/s				
Down Flow Velocity	0.33 \pm 0.025m/s				
HEPA Filter	Two,99.999% efficiency at 0.3um				
Front Window	Motorized.Two-layer laminated toughened glass \geq 5mm				
Noise	EN12469 \leq 58dB/NSF49 \leq 61dB				
UV Lamp	18W*1	18W*1	30W*1	40W*1	40W*1
Fluorescent Lamp	21W*2	21W*2	28W*2	40W*2	40W*1
Consumption	760W	1060W	1300W	1350W	1700W
Illumination	\geq 1000Lux				
Water & Gas Tap	Water Tap*1; Gas Tap*1.(Optional)				
Socket	Two, total consumption:500W				
Display	LCD Display				
Control System	Microprocessor				
Airflow System	70% Air Recirculation, 30% Air Exhaust				
Material	Work Table: 304 stainless steel				
	Main Body: Cold-rolled steel with anti-bacteria powder				
Work Surface Height	750mm				
Caster	Universal wheel with brake and leveling feet				
Power Supply	110/220V+10%. 50/60Hz				
Standard Accessory	Fluorescent lamp,UV lamp*2,Base stand,Remote control,Foot switch,Drain valve, Waterproof sockets				
Optional Accessory	Armrest,Airflow tester,Formanlin fumigation sterilizer,Infrared sterilizer				
Gross Weight	563 lbs	673 lbs	772 lbs	860 lbs	926 lbs
Package Size (W*D*H)(mm)	1250*1050*1870	1450*1050*1870	1640*1040*1860	1940*1050*1960	2090*1050*1860