

Platelet Incubator

Features:

1. Micro-controller controls the internal temperature within $22\pm 2^{\circ}\text{C}$.
2. High accuracy sensor and direct air flow system ensure even temperature distribution.
3. Formed in place CFC free PUF insulation ensures temperature stability and reduced energy consumption.
4. Audio visual alarm for temperature deviation beyond $\pm 2^{\circ}\text{C}$ and power failure.
5. Motion failure alarm when the agitator fails to operate.
6. Automatic agitator and UV lamp stop when door is open.
7. LED display for easy of reading.
8. UV lamp for sterilization.



Model	PLT055	PLT056
Internal Size(W*D*H)	427*470*467 mm	607*535*657 mm
External Size(W*D*H)	540*605*960 mm	720*670*1150 mm
Number of Agitator	Single	
Temperature	$22\pm 2^{\circ}\text{C}$	
Accuracy	$\pm 0.1^{\circ}\text{C}$	
Oscillation	60 ± 5 cycles/minute	
Amplitude	50mm	
Temperature Control	Micro Controller	
Display	LED	
Trays	5	10
Capacity	5-10 plasma bags	10-20 plasma bags
Consumption	101W	180W
Package Size(mm)	620*660*1050 mm	800*725*1240 mm
Net Weight	217 lbs	243 lbs
Power Supply	110V/220V, 50/60Hz	