

BA-D ELECTRONIC DENSITY BALANCE (SPECIFIC GRAVITY)



DETAILS:

BA-D series electronic density (specific gravity) analytic balance is a smart balance specially developed with the adoption of electromagnetic balance sensor and single chip microcomputer of high stability on a basis of BA series electronic balance.

Function

- * Multiple functions: auto-counting, unit conversion (metric carat, gold ounce, etc.), stability, full-scale faring, zero-memory, etc.
- * Built-in automatic calibration guarantees easier, faster, stable and accurate operation of balance.
- * Built-in RS232C output interface guarantees direct connection to computer and printer.
- * Built-in lower weighing hook guarantees convenient use.

MODEL:	PLBA 1091	PLBA 1092	PLBA 1093	PLBA 1094	PLBA 1095
Weighing Range (g)	0-210	0-200	0-100	0-500	0-300
Accuracy (mg)	0.1	0.1	0.1	1	1
Scale size (mm)	φ90	φ90	φ90	φ90	φ90
External Size (mm)	350*215*340	350*215*340	350*215*340	350*215*340	350*215*340
Workspace Height (mm)	240	240	240	240	240
Net Weight (Lbs)	37	37	37	37	37