

Shaking Incubator

Features:

1. Mirror stainless steel chamber.
2. Microprocessor control system, LCD display.
3. Parameter memory function with password lock.
4. Power failure protection and automatic resume.
5. Interlock function: toughed glass cover and shaking motor.
6. Over-temperature protection.



Optional accessories

- a. Desk printer
- b. RS485 interface for software connection

Model	PLT105	PLT106	PLT107
Rotary Speed	20~300rpm		
Rotary Speed Precision	±1 rpm		
Vibration Range	Φ 26mm		
Max Capacity	50ml*12 or 100ml*9	50ml*20 or 100ml*16 or 250ml*12 or 500ml*9	100ml*20 or 250ml*16 or 500ml*12
Standard Configuration	100ml*9	50ml*4, 100ml*4, 250ml*3, 500ml*3	50ml*5, 100ml*5, 250ml*4, 500ml*3
Rocking Plate Size	295*245mm	450*370mm	450*410mm
Timing Range	0~999 hours		
Temperature Range	+5~60℃		4~60℃
Temperature Precision	±0.1℃		
Temperature Fluctuation	±1℃		
Display	LCD Display		
Rocking Plate No.	1 pc(Accessories available)		
Capacity	20L	48L	70L
Consumption	220W	320W	580W
Power Supply	110/220V±10%,50/60Hz		
External Size(W*D*H)	440*410*390mm	600*580*510mm	780*682*560mm
Package Size(W*D*H)	580*530*540mm	740*700*660mm	920*800*610mm
Gross Weight	90.5 lbs	194 lbs	275.5 lbs