

## Drying Oven

### Features:

1. External material: excellent cold-rolled steel with powder coating.
2. Stainless steel inner chamber, round angle structure, adjustable shelves.
3. Double-layer glass observation window.
4. Over-temperature protection.
5. Optional RS485 interface for software connection.



Bench Top Type



Floor Type

|                          |                      |             |             |              |              |               |
|--------------------------|----------------------|-------------|-------------|--------------|--------------|---------------|
| Model                    | PLT150               | PLT151      | PLT152      | PLT153       | PLT154       | PLT155        |
| Capacity                 | 25L                  | 50L         | 70L         | 136L         | 225L         | 429L          |
| Type                     | Bench Top Type       |             | Floor Type  |              |              |               |
| Temperature Range        | RT+10~250°C          |             | RT+10~300°C |              |              |               |
| Temperature Precision    | 0.1°C                |             |             |              |              |               |
| Temperature Fluctuation  | ±0.5°C               |             |             |              |              |               |
| Ambient Temperature      | +5~40°C              |             |             |              |              |               |
| Timing Range             | 1~9999min            |             |             |              |              |               |
| Shelves No.              | 2pcs                 |             |             |              |              |               |
| Consumption              | 500W                 | 750W        | 1240W       | 1740W        | 2500W        | 4500W         |
| Power Supply             | 110/220V±10%,50/60Hz |             |             |              |              | 380V, 50/60Hz |
| Internal Size(W*D*H)(mm) | 300*300*270          | 420*350*350 | 400*350*500 | 450*550*550  | 500*600*750  | 600*550*1300  |
| External Size(W*D*H)(mm) | 500*450*440          | 710*520*520 | 600*550*755 | 600*710*845  | 650*755*1045 | 745*750*1670  |
| Package Size(W*D*H)(mm)  | 730*640*580          | 840*690*670 | 680*680*990 | 730*880*1130 | 780*930*1240 | 880*880*1820  |
| Gross Weight             | 99 lbs               | 110.5 lbs   | 165.5 lbs   | 209.5 lbs    | 242.5 lbs    | 320 lbs       |